


interscience



Anabac[®]

Autoclave deodorant

Anabac[®] Brochure N° 30910-0

GB

Give some air to your lab!

Laboratories are very concerned about hygiene cleanliness and comfort of work. Try **Anabac®**, autoclave deodorant, for a perfumed touch in your lab. **Anabac®** autoclave deodorant bring the discrete final touch, the refinement of cleanliness.



Classic

Well perfumed
All applications

- *mint and eucalyptus fragrance*



Poma

Sweet perfume
Food industry

- *apple fragrance*



Citrus

Sweet perfume
Food industry

- *lemon fragrance*



Floral

Well perfumed
All applications

- *flower fragrance*




Peach

Well perfumed
All applications

- *peach fragrance*




1.




During sample sterilization, bad smell can emerge from the autoclave.

2.



Place one capsule (or 2 if the odor is strong) in the water at the bottom of the autoclave.

3.



After a sterilization cycle, the capsule melts, allowing the perfume to be released.

Advantages

- *Work comfort*
- *Compatible with any autoclave*
- *No effect on the autoclaving*
- *Pleasant perfume*
- *100% biodegradable*

Applications

- For any kind of autoclave (max: 200 L), recommended for food industry.

Specifications

- Pot of 100 capsules
- Pot weight: 160 g
- 5 fragrances:
 - Classic (ref: 320 100)
 - Poma (ref: 320 200)
 - Citrus (ref: 320 300)
 - Floral (ref: 320 400)
 - Peach (ref: 320 500)



Please return this form to your local distributor or to interscience.



Anabac® autoclave deodorant

- ☐ Quotation
- ☐ Free "5 fragrances" sample

I would like to receive more info on:

- ☐ **BagMixer® 400 ml** blender
- ☐ **JumboMix® 3500 ml** blender
- ☐ **BagFilter®** filter bags
- ☐ **BagTools®** accessories
- ☐ **Gravimat®** gravimetric dilutor
- ☐ **easySpiral®** automatic plater
- ☐ **Scan®** colony counters

Name.....

Company

Address

Zip/City..... Country

Phone..... Fax

email.....

web

interscience

Worldwide: 30. Ch Bois Arpents, 78860 St Nom FRANCE
Tel. +33(0)1 34 62 62 61 / Fax. +33(0)1 34 62 43 03
info@interscience.fr / www.interscience.fr

USA/Canada: 199 Weymouth St. ROCKLAND, MA 02370 USA
Tel. +1 781 792 2133 / Fax. +1 781 792 2134
info@intersciencelab.com / www.intersciencelab.com

In addition



Clinet® 101

*Chlorhexidine digluconate
disinfecting wipes*

Ideal for laboratories, hospitals, samplers, home use.



Clinet® 202

*Individually packed chlorhexidine digluconate disinfecting
wipes*

Ideal for laboratories, hospitals, samplers, transports, hand hygiene.



Clinet® 303

Bactericidal and fungicidal spray

Ideal for large surfaces, hard to reach places, angle of walls, rough edges.



Clinet® 404 SaniFood®

*Food contact approved disinfecting
wipes*

For food professionals:
disinfection without rinsing
of surfaces and materials in
contact with food.



Clinet® Anti-algae

*Algicide for laboratories water
baths.*

1 liter bottle with incorporated
doser, Clinet®Anti-algae is used
diluted at 0.1% to inhibit the
growth of algae in water baths.